

HDAC1-specific antibody

Product Information

Catalog No.:	FNab03795
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	histone deacetylase 1
Synonyms:	Histone deacetylase 1 (HD1) Protein deacetylase HDAC1 Protein deacetylase HDAC1 HDAC1 RPD3L1
Observed MW:	60 kDa
Uniprot ID :	Q13547

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:10-1:100
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab03795(HDAC1-specific antibody) at dilution of 1:500