

## HDAC3-specific antibody

### Product Information

Catalog No.:	FNab03799
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	histone deacetylase 3
Synonyms:	Histone deacetylase 3 (HD3) Protein deacetylase HDAC3 Protein deacylase HDAC3 RPD3-2 SMAP45 HDAC3
Observed MW:	50 kDa
Uniprot ID :	O15379

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab03799(HDAC3-specific antibody) at dilution of 1:200