

## HDHC3 antibody

### Product Information

Catalog No.:	FNab03807
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


ppGpp hydrolyzing enzyme involved in starvation response.

### Immunogen information


Immunogen:	HD domain containing 3
Synonyms:	Guanosine-3',5'-bis(diphosphate) 3'-pyrophosphohydrolase MESH1 HD domain-containing protein 3 Metazoan SpoT homolog 1 (MESH1) Penta-phosphate guanosine-3'-pyrophosphohydrolase ((ppGpp)ase) HDHC3 MESH1
Observed MW:	22-25 kDa
Uniprot ID :	Q8N4P3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab03807(HDDC3 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03807(HDDC3 antibody) at dilution of 1:500