

HDGFL3 antibody

Product Information

Catalog No.:	FNab03810
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information


Immunogen:	hepatoma-derived growth factor, related protein 3
Synonyms:	Hepatoma-derived growth factor-related protein 3 (HRP-3) Hepatoma-derived growth factor 2 (HDGF-2) HDGFL3 HDGF2 HDGFRP3
Observed MW:	28 kDa
Uniprot ID :	Q9Y3E1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human uterine cancer using FNab03810(HDGFRP3 antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab03810(HDGFRP3 antibody) at dilution of 1:1000