

## HDHD3 antibody

### Product Information

Catalog No.:	FNab03811
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	haloacid dehalogenase-like hydrolase domain containing 3
Synonyms:	Haloacid dehalogenase-like hydrolase domain-containing protein 3 HDHD3 C9orf158
Observed MW:	28 kDa
Uniprot ID :	Q9BSH5

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab03811(HDHD3 antibody) at dilution of 1:500