

## HDLBP antibody

### Product Information

Catalog No.:	FNab03812
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Appears to play a role in cell sterol metabolism. It may function to protect cells from over-accumulation of cholesterol.

### Immunogen information


Immunogen:	high density lipoprotein binding protein
Synonyms:	Vigilin High density lipoprotein-binding protein (HDL-binding protein) HDLBP HBP VGL
Observed MW:	141 kDa
Uniprot ID :	Q00341

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-HDLBP (IP:FNab03812, 4ug;  
Detection:FNab03812 1:500) with Jurkat cells  
lysate 2400ug.



NIH/3T3 cells were subjected to SDS PAGE  
followed by western blot with FNab03812(HDLBP  
antibody) at dilution of 1:500