

HPN antibody

Product Information

Catalog No.:	FNab03832
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	hepsin
Synonyms:	Serine protease hepsin Transmembrane protease serine 1 Serine protease hepsin non-catalytic chain Serine protease hepsin catalytic chain HPN TMPRSS1
Observed MW:	28 kDa, 45-50 kDa
Uniprot ID :	P05981

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



human liver tissue were subjected to SDS PAGE followed by western blot with FNab03832(Hepsin antibody) at dilution of 1:400