

HERPUD1 antibody

Product Information

Catalog No.:	FNab03838
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	homocysteine-inducible, endoplasmic reticulum stress-inducible, ubiquitin-like domain member 1
Synonyms:	Homocysteine-responsive endoplasmic reticulum-resident ubiquitin-like domain member 1 protein Methyl methanesulfonate (MMF)-inducible fragment protein 1 HERPUD1 HERP KIAA0025 MIF1
Observed MW:	52 kDa
Uniprot ID :	Q15011

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNa03838(HERPUD1 antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNa03838(HERPUD1 antibody) at dilution of 1:600