

HES5 antibody

Product Information

Catalog No.:	FNab03840
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcriptional repressor of genes that require a bHLH protein for their transcription. Plays an important role as neurogenesis negative regulator(By similarity).

Immunogen information


Immunogen:	hairy and enhancer of split 5(Drosophila)
Synonyms:	Transcription factor HES-5 Class B basic helix-loop-helix protein 38 (bHLHb38) Hairy and enhancer of split 5 HES5 BHLHB38
Observed MW:	41 kDa
Uniprot ID :	Q5TA89

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNaB03840(HES5 Antibody) at dilution of 1:50



rat brain tissue were subjected to SDS PAGE followed by western blot with FNaB03840(HES5 Antibody) at dilution of 1:300