

PPIP5K2 antibody

Product Information

Catalog No.:	FNab03873
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	histidine acid phosphatase domain containing 1
Synonyms:	Inositol hexakisphosphate and diphosphoinositol-pentakisphosphate kinase 2 Diphosphoinositol pentakisphosphate kinase 2 Histidine acid phosphatase domain-containing protein 1 InsP6 and PP-IP5 kinase 2 VIP1 homolog 2 (hsVIP2) PPIP5K2 HISPPD1 KIAA0433 VIP2
Observed MW:	138-140 kDa
Uniprot ID :	O43314

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab03873(HISPPD1 Antibody) at dilution of 1:1000