

H4C1 antibody

Product Information

Catalog No.:	FNab03878
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	histone cluster 1, H4f
Synonyms:	Histone H4 H4C1 H4/A H4FA HIST1H4A
Observed MW:	11 kDa
Uniprot ID :	P62805

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



human spleen tissue were subjected to SDS PAGE followed by western blot with FNab03878(HIST1H4F antibody) at dilution of 1:300