

HLA-B antibody

Product Information

Catalog No.:	FNab03897
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	major histocompatibility complex, class I, B
Synonyms:	HLA class I histocompatibility antigen, B alpha chain Human leukocyte antigen B (HLA-B) HLA-B HLAB
Observed MW:	40-45 kDa
Uniprot ID :	P18464

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:100-1:4000; IP: 1:1000-1:4000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	




Immunohistochemistry of paraffin-embedded human spleen using FNab03897(HLA class I(HLA-B) antibody) at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using FNab03897(HLA-B antibody) at dilution of 1:50



IP Result of anti-HLA class I (HLA-B)
(IP:FNab03897, 4ug; Detection:FNab03897
1:2000) with HepG2 cells lysate 4000ug.



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03897(HLA class I(HLA-B) antibody) at dilution of 1:800