

## HLA-F antibody

### Product Information

Catalog No.:	FNab03911
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the presentation of foreign antigens to the immune system.

### Immunogen information


Immunogen:	major histocompatibility complex, class I, F
Synonyms:	HLA class I histocompatibility antigen, alpha chain F CDA12 HLA F antigen Leukocyte antigen F MHC class I antigen F HLA-F HLA-5.4 HLAF
Observed MW:	45 kDa
Uniprot ID :	P30511

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue using FNab03911(HLA-F antibody) at dilution of 1:100



A431 cells were subjected to SDS PAGE followed by western blot with FNab03911(HLA-F antibody) at dilution of 1:500