

HMGA2 antibody

Product Information

Catalog No.:	FNab03923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Functions as a transcriptional regulator. Functions in cell cycle regulation through CCNA2. Plays an important role in chromosome condensation during the meiotic G2/M transition of spermatocytes.

Immunogen information

Immunogen:	high mobility group AT-hook 2
Synonyms:	High mobility group protein HMGI-C High mobility group AT-hook protein 2 HMGA2 HMGIC
Observed MW:	18-20 kDa
Uniprot ID :	P52926

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab03923(HMGA2 Antibody) at dilution of 1:1000