

## HMGN2 antibody

### Product Information

Catalog No.:	FNab03932
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds to the inner side of the nucleosomal DNA thus altering the interaction between the DNA and the histone octamer. May be involved in the process which maintains transcribable genes in a unique chromatin conformation(By similarity).

### Immunogen information


Immunogen:	high-mobility group nucleosomal binding domain 2
Synonyms:	Non-histone chromosomal protein HMG-17 High mobility group nucleosome-binding domain-containing protein 2 HMGN2 HMG17
Observed MW:	18-20 kDa
Uniprot ID :	P05204

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab03932(HMGN2 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03932(HMGN2 antibody) at dilution of 1:600