

## HMGN4 antibody

### Product Information

Catalog No.:	FNab03933
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	high mobility group nucleosomal binding domain 4
Synonyms:	High mobility group nucleosome-binding domain-containing protein 4 Non-histone chromosomal protein HMG-17-like 3 (Non-histone chromosomal protein) HMGN4 HMG17L3 NHC
Observed MW:	15 kDa
Uniprot ID :	O00479

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03933(HMGN4 antibody) at dilution of 1:400