

## HMGXB4 antibody

### Product Information

Catalog No.:	FNab03934
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Negatively regulates Wnt/beta-catenin signaling during development.

### Immunogen information

Immunogen:	HMG box domain containing 4
Synonyms:	HMG domain-containing protein 4 HMG box-containing protein 4 High mobility group protein 2-like 1 Protein HMGBCG HMGXB4 HMG2L1 HMGBCG
Observed MW:	66 kDa
Uniprot ID :	Q9UGU5

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:1000-1:4000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03934( HMGXB4 Antibody) at dilution of 1:2000