

HMMR-Specific antibody

Product Information

Catalog No.:	FNab03936
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	hyaluronan-mediated motility receptor(RHAMM)
Synonyms:	Hyaluronan mediated motility receptor Intracellular hyaluronic acid-binding protein Receptor for hyaluronan-mediated motility HMMR IHABP RHAMM
Observed MW:	70 kDa
Uniprot ID :	O75330

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab03936(HMMR-Specific Antibody) at dilution of 1:200



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03936(HMMR-Specific antibody) at dilution of 1:800