

## HNF1A antibody

### Product Information

Catalog No.:	FNab03941
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcriptional activator that regulates the tissue specific expression of multiple genes, especially in pancreatic islet cells and in liver. Required for the expression of several liver specific genes. Binds to the inverted palindrome 5'-GTTAATNATTAAC-3'.

### Immunogen information

Immunogen:	HNF1 homeobox A
Synonyms:	Hepatocyte nuclear factor 1-alpha (HNF-1-alpha, HNF-1A) Liver-specific transcription factor LF-B1 (LFB1) Transcription factor 1 (TCF-1) HNF1A TCF1
Observed MW:	67 kDa
Uniprot ID :	P20823

### Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



rat liver tissue were subjected to SDS PAGE followed by western blot with FNab03941(HNF1A antibody) at dilution of 1:1000