

HNF1B antibody

Product Information

Catalog No.:	FNab03942
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcription factor, probably binds to the inverted palindrome 5'-GTTAATNATTAAC-3'.

Immunogen information


Immunogen:	HNF1 homeobox B
Synonyms:	Hepatocyte nuclear factor 1-beta (HNF-1-beta, HNF-1B) Homeoprotein LFB3 Transcription factor 2 (TCF-2) Variant hepatic nuclear factor 1 (vHNF1) HNF1B TCF2
Observed MW:	58-60 kDa
Uniprot ID :	P35680

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab03942(HNF1B Antibody) at dilution of 1:400



L02 cells were subjected to SDS PAGE followed by western blot with FNab03942(HNF1B antibody) at dilution of 1:800