

## HNF4G antibody

### Product Information

Catalog No.:	FNab03943
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcription factor. Has a lower transcription activation potential than HNF4-alpha.

### Immunogen information


Immunogen:	hepatocyte nuclear factor 4, gamma
Synonyms:	Hepatocyte nuclear factor 4-gamma (HNF-4-gamma) Nuclear receptor subfamily 2 group A member 2 HNF4G NR2A2
Observed MW:	46 kDa
Uniprot ID :	Q14541

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using FNab03943(HNF4G Antibody) at dilution of 1:50



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab03943(HNF4G Antibody) at dilution of 1:600