

HNRNPA3 antibody

Product Information

Catalog No.:	FNab03948
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in cytoplasmic trafficking of RNA. Binds to the cis-acting response element, A2RE. May be involved in pre-mRNA splicing.

Immunogen information


Immunogen:	heterogeneous nuclear ribonucleoprotein A3
Synonyms:	Heterogeneous nuclear ribonucleoprotein A3 (hnRNP A3) HNRNPA3 HNRPA3
Observed MW:	40 kDa
Uniprot ID :	P51991

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunofluorescent analysis of HeLa cells using FNab03948 (HNRNPA3 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03948(HNRNPA3 Antibody) at dilution of 1:1000