

## HNRNPH3 antibody

### Product Information

Catalog No.:	FNab03954
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the splicing process and participates in early heat shock-induced splicing arrest. Due to their great structural variations the different isoforms may possess different functions in the splicing reaction.

### Immunogen information

Immunogen:	heterogeneous nuclear ribonucleoprotein H3(2H9)
Synonyms:	Heterogeneous nuclear ribonucleoprotein H3 (hnRNP H3) Heterogeneous nuclear ribonucleoprotein 2H9 (hnRNP 2H9) HNRNPH3 HNRPH3
Observed MW:	37-40 kDa
Uniprot ID :	P31942

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab03954( HNRNPH3 Antibody) at dilution of 1:300