

HNRNPR antibody

Product Information

Catalog No.:	FNab03957
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


This gene encodes an RNA-binding protein that is a member of the spliceosome C complex, which functions in pre-mRNA processing and transport. The encoded protein also promotes transcription at the c-fos gene. Alternative splicing results in multiple transcript variants. There are pseudogenes for this gene on chromosomes 4, 11, and 10.

Immunogen information


Immunogen:	heterogeneous nuclear ribonucleoprotein R
Synonyms:	Heterogeneous nuclear ribonucleoprotein R (hnRNP R) HNRNPR HNRPR
Observed MW:	70 kDa
Uniprot ID :	O43390

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:100 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded rat brain using FNab03957(HNRNPR antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03957(HNRNPR antibody) at dilution of 1:1000