

RNF31 antibody

Product Information

Catalog No.:	FNab03962
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ring finger protein 31
Synonyms:	E3 ubiquitin-protein ligase RNF31 HOIL-1-interacting protein (HOIP) RING finger protein 31 RING-type E3 ubiquitin transferase RNF31 Zinc in-between-RING-finger ubiquitin-associated domain protein RNF31 ZIBRA
Observed MW:	100 kDa
Uniprot ID :	Q96EP0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab03962(HOIP antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03962(HOIP antibody at dilution of 1:200