

## HOXB1 antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab03977  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Acts on the anterior body structures.

### Immunogen information

|              |                |
|--------------|----------------|
| Immunogen:   | homeobox B1    |
| Synonyms:    | HOXB1 antibody |
| Observed MW: | 38-44 kDa      |
| Uniprot ID : | P14653         |

### Application

|                          |              |
|--------------------------|--------------|
| Reactivity:              | Human, Mouse |
| Tested Application:      | ELISA, WB    |
| Recommended dilution:WB: | 1:500-1:2000 |
| Image:                   |              |



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03977( HOXB1 Antibody) at dilution of 1:300