

HOXB7 antibody

Product Information

Catalog No.:	FNab03978
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Immunogen information

Immunogen:	homeobox B7
Synonyms:	Homeobox protein Hox-B7 Homeobox protein HHO.C1 Homeobox protein Hox-2C HOXB7 HOX2C
Observed MW:	25-30 kDa
Uniprot ID :	P09629

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



SKOV-3 cells were subjected to SDS PAGE followed by western blot with FNab03978(HOXB7 antibody) at dilution of 1:1500