

## HOXD4 antibody

### Product Information

Catalog No.:	FNab03988
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

### Immunogen information

Immunogen:	homeobox D4
Synonyms:	Homeobox protein Hox-D4 Homeobox protein HHO.C13 Homeobox protein Hox-4B Homeobox protein Hox-5.1 HOXD4 HOX4B
Observed MW:	28 kDa
Uniprot ID :	P09016

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab03988(HOXD4 antibody) at dilution of 1:1000