

HOXD8 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab03989 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Immunogen information

| | |
|--------------|--|
| Immunogen: | homeobox D8 |
| Synonyms: | Homeobox protein Hox-D8 Homeobox protein Hox-4E Homeobox protein Hox-5.4 HOXD8 HOX4E |
| Observed MW: | 32 kDa |
| Uniprot ID : | P13378 |

Application

| | |
|---------------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution: WB: | 1:500-1:2000 |
| Image: | |



A375 cells were subjected to SDS PAGE followed by western blot with FNab03989(HOXD8 Antibody) at dilution of 1:300