

PLAAT5 antibody

Product Information

Catalog No.:	FNab04007
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


catalyzes N-acylation of phosphatidylethanolamine(PE) to generate N-Acylphosphatidylethanolamine(NAPE) a precursor of bioactive N-acylethanolamines, including the endocannabinoid anandamide.

Immunogen information

Immunogen:	HRAS-like suppressor family, member 5
Synonyms:	Phospholipase A and acyltransferase 5 Ca(2+)-independent N-acyltransferase (iNAT) H-rev107-like protein 5 HRAS-like suppressor 5 (HRSL5) PLAAT5 HRASLS5 HRLP5
Observed MW:	None
Uniprot ID :	Q96KN8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNa04007(HRASLS5 antibody) at dilution of 1:100