

HRAS-Specific antibody

Product Information

Catalog No.:	FNab04008
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	v-Ha-ras Harvey rat sarcoma viral oncogene homolog
Synonyms:	GTPase HRas H-Ras-1 Ha-Ras Transforming protein p21 c-H-ras p21ras GTPase HRas, N-terminally processed HRAS HRAS1
Observed MW:	21 kDa
Uniprot ID :	P01112

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, IHC, WB, FC
Recommended dilution:	WB: 1:200-1:1000; IP:1:1000-1:6000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human kidney slide using FNab04008(HRAS-Specific Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab04008(HRAS-Specific Antibody) at dilution of 1:300