

HSF5 antibody

Product Information

Catalog No.:	FNab04040
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May act as a transcriptional factor.

Immunogen information

Immunogen:	heat shock transcription factor family member 5
Synonyms:	Heat shock factor protein 5 (HSF 5) Heat shock transcription factor 5 (HSTF 5) HSF5 HSTF5
Observed MW:	65 kDa
Uniprot ID :	Q4G112

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab04040(HSF5 antibody) at dilution of 1:2000