

HSPB6 antibody

Product Information

Catalog No.:	FNab04043
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	heat shock protein, alpha-crystallin-related, B6
Synonyms:	Heat shock protein beta-6 (HspB6) Heat shock 20 kDa-like protein p20 HSPB6
Observed MW:	18-20 kDa
Uniprot ID :	O14558

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab04043(HSP20 antibody) at dilution of 1:4000