

## HSPA1A antibody

### Product Information

Catalog No.:	FNab04049
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	heat shock 70kDa protein 1B
Synonyms:	Heat shock 70 kDa protein 1A Heat shock 70 kDa protein 1 (HSP70-1, HSP70.1) HSPA1A HSP72 HSPA1 HSX70
Observed MW:	66-70 kDa
Uniprot ID :	P0DMV8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB,IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded rat brain tissue slide using FNab04048(HSP70 Antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04049( HSP70 Antibody) at dilution of 1:1000