

HSPG2-Specific antibody

Product Information

Catalog No.:	FNab04063
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	heparan sulfate proteoglycan 2
Synonyms:	Basement membrane-specific heparan sulfate proteoglycan core protein (HSPG) Perlecan (PLC) Endorepellin LG3 peptide HSPG2
Observed MW:	469 kDa
Uniprot ID :	P98160

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab04063(HSPG2-Specific antibody) at dilution of 1:300