

TRMT2A antibody

Product Information

Catalog No.:	FNab04066
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	TRM2 tRNA methyltransferase 2 homolog A
Synonyms:	tRNA (uracil-5-)-methyltransferase homolog A mRNA (uracil-5-)-methyltransferase TRMT2A TRMT2A HTF9C
Observed MW:	42-44, 68-75 kDa
Uniprot ID :	Q8IZ69

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab04066(TRMT2A antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab04066(TRMT2A Antibody) at dilution of 1:600