

## ELAVL3 antibody

### Product Information

Catalog No.:	FNab04072
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	ELAV(embryonic lethal, abnormal vision, Drosophila)-like 3(Hu antigen C)
Synonyms:	ELAV-like protein 3 Hu-antigen C (HuC) Paraneoplastic cerebellar degeneration-associated antigen Paraneoplastic limbic encephalitis antigen 21 ELAVL3 HUC PLE21
Observed MW:	36-38 kDa
Uniprot ID :	Q14576

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab04072(ELAVL3 Antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab04072(HuC antibody) at dilution of 1:200