

HuA/B/C/D antibody

Product Information

Catalog No.:	FNab04073
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ELAV(embryonic lethal, abnormal vision, Drosophila)-like 4(Hu antigen D)
Synonyms:	ELAV-like protein 4 Hu-antigen D (HuD) Paraneoplastic encephalomyelitis antigen HuD ELAVL4 HUD PNEM
Observed MW:	40 kDa
Uniprot ID :	P26378

Application

Reactivity:	Human, Mouse, Rat, Zebrafish
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:1000-1:4000; IP: 1:200-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab04073(HuA/B/C/D antibody at dilution of 1:200



human cerebellum tissue were subjected to SDS PAGE followed by western blot with FNab04073(HuA/B/C/D antibody) at dilution of 1:800