

HuD antibody

Product Information

Catalog No.:	FNab04074
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ELAV(embryonic lethal, abnormal vision, Drosophila)-like 4(Hu antigen D)
Synonyms:	ELAV-like protein 4 Hu-antigen D (HuD) Paraneoplastic encephalomyelitis antigen HuD ELAVL4 HUD PNEM
Observed MW:	40 kDa
Uniprot ID :	P26378

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:1000-1:4000; IF: 1:10-1:100
Image:	



Immunofluorescent analysis of SH-SY5Y cells using FNa04074 (HuD antibody) at dilution of 1:25



human cerebellum tissue were subjected to SDS PAGE followed by western blot with FNa04074(HuD antibody) at dilution of 1:800