

MT-RNR2 antibody

Product Information

Catalog No.:	FNab04078
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	l-rRNA
Synonyms:	Humanin Humanin mitochondrial (HNM) MT-RNR2 HN
Observed MW:	2 kDa
Uniprot ID :	Q8IVG9

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab04078(humanin Antibody) at dilution of 1:50