

TNFRSF14 antibody

Product Information

Catalog No.:	FNab04086
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	tumor necrosis factor receptor superfamily, member 14(herpesvirus entry mediator)
Synonyms:	Tumor necrosis factor receptor superfamily member 14 Herpes virus entry mediator A (Herpesvirus entry mediator A, HveA) Tumor necrosis factor receptor-like 2 (TR2) TNFRSF14 HVEA HVEM
Observed MW:	None
Uniprot ID :	Q92956

Application

Reactivity:	Human
Tested Application:	ELISA