

IFFO1 antibody

Product Information

Catalog No.:	FNab04132
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	intermediate filament family orphan 1
Synonyms:	Non-homologous end joining factor IFFO1 (NHEJ factor IFFO1) Intermediate filament family orphan 1 Tumor antigen HOM-TES-103 IFFO1 IFFO
Observed MW:	75-85 kDa
Uniprot ID :	Q0D2I5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNaB04132 (IFFO1 antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNaB04132 (IFFO1 Antibody) at dilution of 1:1000