

IFIT1B antibody

Product Information

Catalog No.:	FNab04136
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	interferon-induced protein with tetratricopeptide repeats 1-like
Synonyms:	Interferon-induced protein with tetratricopeptide repeats 1B Interferon-induced protein with tetratricopeptide repeats 1-like protein IFIT1B IFIT1L
Observed MW:	60 kDa
Uniprot ID :	Q5T764

Application

Reactivity:	Human, Rat, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200 - 1:1000; IHC: 1:20 - 1:200; IF: 1:10 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab04136(IFIT1L Antibody) at dilution of 1:50



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with FNab04136(IFIT1L Antibody) at dilution of 1:300