

IFT25 antibody

Product Information

Catalog No.:	FNab04156
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	heat shock protein family B(small), member 11
Synonyms:	Intraflagellar transport protein 25 homolog Heat shock protein beta-11 (Hspb11) Placental protein 25 (PP25) IFT25 C1orf41 HSPB11
Observed MW:	25 kDa
Uniprot ID :	Q9Y547

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500 - 1:2400;IF: 1:20 - 1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04156 (IFT25 antibody) at dilution of 1:500