

IGF1B-Specific antibody

Product Information

Catalog No.:	FNab04169
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	insulin-like growth factor 1(somatomedin C)
Synonyms:	Insulin-like growth factor I (IGF-I) Mechano growth factor (MGF) Somatomedin-C IGF1 IBP1
Observed MW:	22 kDa
Uniprot ID :	P05019

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:20 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab04169(IGF1B-Specific antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab04169(IGF1B-Specific antibody) at dilution of 1:500