

MIIP antibody

Product Information

Catalog No.:	FNab04195
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	invasion inhibitory protein 45
Synonyms:	Migration and invasion-inhibitory protein IGFBP2-binding protein Invasion-inhibitory protein 45 (Iip45) MIIP IIP45
Observed MW:	45 kDa
Uniprot ID :	Q5JXC2

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab04195(IIP45 antibody) at dilution of 1:200



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab04195(IIP45 Antibody) at dilution of 1:600