

NFKBIA antibody

Product Information

Catalog No.:	FNab04198
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, alpha
Synonyms:	NF-kappa-B inhibitor alpha I-kappa-B-alpha (IκB-alpha, IκappaBalpha) Major histocompatibility complex enhancer-binding protein MAD3 NFKBIA IKBA MAD3 NFKBI
Observed MW:	39 kDa
Uniprot ID :	P25963

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:50-1:100
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab04198(IKBA Antibody) at dilution of 1:100



A549 cells were subjected to SDS PAGE followed by western blot with FNab04198(IKBA antibody) at dilution of 1:1000