

IKBIP antibody

Product Information

Catalog No.:	FNab04203
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	IKK interacting protein
Synonyms:	Inhibitor of nuclear factor kappa-B kinase-interacting protein (I kappa-B kinase-interacting protein, IKBKB-interacting protein, IKK-interacting protein) IKBIP IKIP
Observed MW:	40 kDa
Uniprot ID :	Q70UQ0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab04203(IKIP antibody) at dilution of 1:500