

IKZF3 antibody

Product Information

Catalog No.:	FNab04207
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcription factor that plays an important role in the regulation of lymphocyte differentiation. Plays an essential role in regulation of B-cell differentiation, proliferation and maturation to an effector state. Involved in regulating BCL2 expression and controlling apoptosis in T-cells in an IL2-dependent manner.

Immunogen information


Immunogen:	IKAROS family zinc finger 3(Aiolos)
Synonyms:	Zinc finger protein Aiolos Ikaros family zinc finger protein 3 IKZF3 ZNFN1A3
Observed MW:	58-70 kDa
Uniprot ID :	Q9UKT9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human spleen using FNab04207(IKZF3 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab04207(IKZF3 antibody) at dilution of 1:300