

IL12A antibody

Product Information

Catalog No.:	FNab04218
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3C7
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated Killer cells, and stimulate the production of IFN-gamma by resting PBMC.

Immunogen information


Immunogen:	interleukin 12A(natural killer cell stimulatory factor 1, cytotoxic lymphocyte maturation factor 1, p35)
Synonyms:	Interleukin-12 subunit alpha (IL-12A) Cytotoxic lymphocyte maturation factor 35 kDa subunit (CLMF p35) IL-12 subunit p35 NK cell stimulatory factor chain 1 (NKSF1) IL12A NKSF1
Observed MW:	25 kDa
Uniprot ID :	P29459

Application

Reactivity:	Human
Tested Application:	ELISA, WB,IHC
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab04218(IL12A Antibody) at dilution of 1:200



Recombinant protein were subjected to SDS PAGE followed by western blot with FNab04218(IL12A Antibody) at dilution of 1:2000